



Zeiss Stemi 305

mit PL 10x 23 Br. foc. Okular | with PL 10x 23 Br. foc. eyepiece

Vorsatzoptik <i>Front optic</i>	Arbeitsabstand <i>Free working distance</i>		Objektfeld <i>Object field</i>
0,5x	185 mm	4x – 20x	57,5 – 11,5 mm
0,75x	128 mm	6x – 30x	38,3 – 7,7 mm
1x ohne Vorsatzoptik <i>without front optic</i>	110 mm	8x – 40x	28,8 – 5,8 mm
1,5x	56 mm	12x – 60x	19,2 – 3,8 mm
2x	43 mm	16x – 80x	14,4 – 2,9 mm

Zeiss Stemi 305

mit PL 16x 14 Br. foc. Okular | with PL 16x 14 Br. foc. eyepiece

Vorsatzoptik <i>Front optic</i>	Arbeitsabstand <i>Free working distance</i>		Objektfeld <i>Object field</i>
0,5x	185 mm	6,4x – 32x	35 – 7 mm
0,75x	128 mm	9,6x – 48x	23,3 – 4,7 mm
1x ohne Vorsatzoptik <i>without front optic</i>	110 mm	12,8x – 64x	17,5 – 3,5 mm
1,5x	56 mm	19,2x – 96x	11,7 – 2,3 mm
2x	43 mm	25,6x – 128x	8,8 – 1,8 mm